

NOW AVAILABLE  
IN WHITE!

# Bed-r'Nest®

PORTION-CONTROLLED ENRICHMENT  
AND NESTING MATERIAL



## A NATURE-LIKE NEST FROM DAY ONE

### PORTION CONTROLLED AND EASY TO DISPENSE

Available in 4 gram and 8 gram sizes, Bed-r'Nest is a portion-controlled, easy-to-dispense enrichment and nesting material.

### ENHANCED CAGE ENVIRONMENT

Bed-r'Nest is designed to resemble the thin strands of grass and other materials mice instinctively use to build nests in nature. With no outer shell or other material to work through, mice have immediate access to enrichment and nesting material.

Bed-r'Nest also aids mice in thermoregulation, promoting less stress in the animals and, ultimately, better research.

### LABOR SAVINGS

As a portion-controlled, easy-to-dispense enrichment and nesting solution, Bed-r'Nest eliminates the difficulties associated with the handling and portioning of bulk nesting material.

With no pulling, tearing or prepping required, Bed-r'Nest can be removed from packaging and immediately dispensed into cages.

### CONSISTENT, PREDICTABLE COST PER CAGE

Bed-r'Nest provides a consistent amount of enrichment and nesting material, so you can count on a predictable cost per cage, every time.

### NOW AVAILABLE IN WHITE

White paper aids in health checks and seeing pups in the nest, and is perfect for toxicology studies. No chlorine bleach is used in the manufacturing of this product.

### PRODUCT SPECIFICATIONS

Bed-r'Nest is available with either 4 grams or 8 grams of enrichment and nesting material. Bed-r'Nest is available in both kraft and white paper, and in cases with sleeves or in bulk. Bed-r'Nest is available irradiated only in sleeved SKUs. Bulk SKUs are only available non-irradiated. For complete product listings, see the reverse of this sheet.

## Bed-r'Nest®

Consistent, portion-controlled enrichment and nesting material ✓

Arrives ready-to-use out of the box; no tearing or separating required ✓

Allows for easy-to-perform health checks ✓

Aids in thermoregulation ✓

Enables the building of a dome-shaped nest ✓

Minimizes fiber residue in cage environment ✓

Available irradiated (sleeved SKUs only) ✓

Available in 4 gram and 8 gram sizes ✓

Available in kraft and white ✓

Available in sleeves and in bulk ✓

# PROVIDING A COMPLETE LINE OF BEDDING, NESTING AND ENRICHMENT PRODUCTS

	PRODUCT	DESCRIPTION	PRODUCT CODE	CAGE FILLS	PACKAGING	PALLET STACK
BEDDING	<b>Bed-o'Cobs® ¼"</b>	1/4" Cob*	4B	115	1.25 cu ft/35 L bag	70
	<b>Bed-o' Cobs® ⅛"</b>	1/8" Cob	8B	115	1.25 cu ft/35 L bag	60
	<b>Bed-o' Cobs® Combo</b>	1/8" & 1/4 Cob	CB	115	1.25 cu ft/35 L bag	70
	<b>Pure-o' Cel®</b>	Paper Squares	PB	115	1.25 cu ft/35 L bag	80
	<b>XTRAFine™</b>	Finely Ground Extra Absorbent Cob	XFB	115	1.25 cu ft/35 L bag	60
NESTING AND ENRICHMENT	<b>Bed-r'Nest®</b> Portion Controlled Nesting	4 gram (kraft)	BRN4SR BRN4BP	1600 3200	1600/cs (40-ct sleeves) 3200/cs (without sleeves)	40 16
		4 gram (white)	BRN4WSR BRN4WBP	1600 3200	1600/cs (40-ct sleeves) 3200/cs (without sleeves)	40 16
		8 gram (kraft)	BRN8SR BRN8BP	800 2200	800/cs (20-ct sleeves) 2200/cs (without sleeves)	40 16
		8 gram (white)	BRN8WSR BRN8WBP	800 2200	800/cs (20-ct sleeves) 2200/cs (without sleeves)	40 16
		<b>Enrich-n'Nest®</b>	Paper Rolls	EB	1600	1.25 cu ft/35 L bag
	<b>Crink-l'Nest™</b>	Crinkled Paper (kraft or white)	CNK CNW	1600-2000	22 lb case	15
	BEDDING AND ENRICHMENT	<b>Bed-o' Cobs® / Pure-o' Cel®</b>	1/4" Cob & Paper Squares	4PB	115	1.25 cu ft/35 L bag
<b>Enrich-o'Cobs®</b>		1/4" Cob & Paper Rolls	4EB	115	1.25 cu ft/35 L bag	80
<b>Enrich-n' Pure®</b>		Paper Rolls & Squares	EPB	115	1.25 cu ft/35 L bag	80
SPECIALTY	<b>LabFresh®</b>	Extra Absorbent Litter	LF20	-	20 lb bag	77
	<b>Pel-o'Cobs®</b>	Pelleted Cob Bedding	PELB	-	1.25 cu ft/35 L bag	50
	<b>Alfa-a'Cobs®</b>	Pelleted Deodorizing Cob Bedding	ALFB	-	1.25 cu ft/35 L bag	50

\*Also available qualified

For more information visit [AndersonsLabBedding.com](http://AndersonsLabBedding.com)

©2017 The Andersons, Inc. All rights reserved. Alf-a'Cobs, The Andersons logo, Bed-o'Cobs, Bed-r'Nest, Enrich-n'Nest, Enrich-n'Pure, Enrich-o'Cobs, LabFresh, Pel-o'Cobs, and Pure-o'Cel are registered trademarks of The Andersons, Inc. Crink-l'Nest, and XTRAFine are trademarks of The Andersons, Inc. H17

