



## **BSE - CB Series Automatic Bedding Filling Machine**

# Automatic bedding filler

The BSE-CB series automatic bedding filling machine is specially used for the automatic filling of animal cage bedding, Including bedding storage, automatic filling, automatic dust removal. To solve the problems of high labor intensity, heavy dust pollution, uneven filling and so on. It is suitable for all kinds of freefall bedding, such as corn cob, etc.



- **Supermatic**

Automatic filling, automatic detection cage, automatic quantitative filling, automatic dust removal, no manual intervention.

- **High efficiency**

15-25 / min cage filling speed. Meet most animal lab needs.

- **Induction blanking**

Automatic detection of cage, automatic start filling process.

- **Ensure the welfare of operators**

The operating area is a negative pressure environment, which can automatically absorb cushion dust during operation to protect the safety of personnel.

- **Easy to operate**

Touch screen operation, inside a variety of cage bedding filling procedures, one-button operation, easy to use.

- **Wide application**

Preset multiple bedding filling procedures, can adapt to various specifications of cage box pad material adding



name	Model	Overall dimensions ( L*W*H )	Power supply	Fill efficiency
Automatic bedding filler	BSE-CB1100	1100 x 740 x 1780mm	380V/50HZ/2kw	49L/min
filling speed	deviation	cleanliness	Air volume	noise
15-25个/min	<20g	Cleanliness level 7	1060m³/h	60dB

## Developed Assembled & Manufactured by

ATNT Laboratories in technical collaboration with

Shinva Medical P.R. China

ISPA engineering Pvt Ltd India

Wagle Industrial Estate, Mumbai

Thane -400604 Maharashtra India

## ATNT Laboratories

Unit # 812, Excellencia Lodha Supremus 1,

Rd No.22, Wagle Estate, Mumbai, Thane - 400604, India.

Mob. No. +91 - 9892520959

Tel. +91-22-25830958 / + 91-22-25830959

Email ID : info@atntlabs.com, ashutosh@atntlabs.com

Web : www.atntlabs.com